



## Hanna Lab App

Available on iOS and Android

The first app that turns a smart phone or tablet into a full-featured pH meter.


The Hanna Lab App turns a compatible smart phone or tablet into a full-featured pH meter when used with HALO®. Functions include calibration, measurement, data logging, graphing, and data sharing. Measurement and logging of pH and temperature at one second intervals start as soon as the probe is connected. Measurements can be displayed alone, with tabulated data, or as a graph. The graph can be panned and zoomed with pinch-to-zoom technology.





## Data Logging

pH	mV	T (°C)	Date	Note	Status
7.212	-11.8	22.9	1/10/17, 2:01:01 PM		
7.211	-11.9	23.0	1/10/17, 2:01:02 PM		
7.209	-12.2	23.0	1/10/17, 2:01:03 PM		
7.207	-12.1	23.0	1/10/17, 2:01:04 PM		
7.207	-11.7	22.9	1/10/17, 2:01:05 PM		
7.200	-11.7	23.0	1/10/17, 2:01:06 PM		
7.202	-11.9	23.0	1/10/17, 2:01:07 PM		
7.200	-11.7	23.0	1/10/17, 2:01:08 PM		
7.202	-11.8	23.0	1/10/17, 2:01:09 PM		
7.203	-11.8	23.0	1/10/17, 2:01:10 PM		
7.203	-11.8	23.0	1/10/17, 2:01:11 PM		
7.203	-11.7	23.0	1/10/17, 2:01:12 PM		
7.200	-12.3	23.0	1/10/17, 2:01:13 PM		
7.202	-12.2	23.0	1/10/17, 2:01:14 PM		
7.203	-11.6	23.0	1/10/17, 2:01:15 PM		
7.203	-11.8	23.0	1/10/17, 2:01:16 PM		
7.203	-11.8	23.0	1/10/17, 2:01:17 PM		
7.203	-11.8	23.0	1/10/17, 2:01:18 PM		
7.200	-11.7	23.0	1/10/17, 2:01:19 PM		
7.200	-11.8	23.0	1/10/17, 2:01:20 PM		
7.199	-11.7	23.0	1/10/17, 2:01:21 PM		
7.200	-11.8	23.0	1/10/17, 2:01:22 PM		
7.200	-11.7	23.0	1/10/17, 2:01:23 PM		
7.200	-11.7	23.0	1/10/17, 2:01:24 PM		



## Data-logging

Data is automatically saved every hour. There are four ways to save and share data: All data since last auto save, Annotations only, All data within a timed interval, and Annotations within a timed interval.

## Export Data

Share data via email in PDF or CSV format.

## Custom Annotations

Saved data points may be annotated with measurement specific information.



## GLP (Good Laboratory Practice)



### Basic GLP

Displays date and time of current calibration along with probe offset and average slope. For tablet displays, basic GLP can be also displayed in table and graph views.



### Full GLP

Displays date and time of current calibration, probe offset, and average slope along with calibrated buffers, mV values, temperature and slopes between each buffer. For tablet displays, full GLP can be also displayed in table and graph views.

### Hanna Lab App Specifications\*

Range**	-2.000 to 16.000 pH; ±800 mV; -20.0 to 120.0°C (-4.0 to 248.0°F)
Resolution	0.1; 0.01; 0.001 pH; 1; 0.1 mV; 0.1°C (0.1°F)
Accuracy (@25°C/77°F)	±0.005 pH; ±0.3 mV; ±0.5°C (±1.0°F)
Calibration Points	up to five-point calibration with seven standard buffers (1.68, 3.00 (HI10482 only) or 4.01, 6.86, 7.01, 9.18, 10.01, 12.45 pH)
Temperature Compensation**	automatic from -5.0 to 100.0 °C – 23.0 to 212.0 °F
Compatibility/System Requirements	see <a href="http://www.hannainst.com">www.hannainst.com</a> for latest compatibility requirements

### Download Information



\* \*\* Limits will be reduced to actual probe/sensor limits.


Apple, the Apple logo, iPhone and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. Android, Google Play and the Google Play logo are trademarks of Google Inc.

## Calibration



### Clear and Concise Calibration Screens


The Hanna Lab App allows for calibration of up to five points. The buffer value is automatically detected and temperature corrected to 25.0°C during calibration.



### Calibration Reminder

Alerts users when HALO needs calibration.

## Additional Features



### Measurement Alerts

Readings that exceed user-defined alarm thresholds are highlighted in yellow on the measurement screen, graph, and table. Readings that exceed the probe specifications are highlighted in red.



### Settings

Tap the gear icon in the top right corner of the measurement screen to access the Settings menu.



### Help and Tutorials:

The Hanna Lab App also features demo probe mode, general app information, general HALO information, pH tutorial, maintenance tutorial, and contact information.