



## Ideal for meat testing

The FC232D pH electrode is recommended for use with the Hanna HI-99163N portable pH meter and is suitable for pH measurement of meat products. This probe has a 1m cable and DIN connector.

The electrode can be tipped with a stainless steel knife blade for meat processing.

Please note: This electrode is not supplied with the stainless steel knife blade for meat processing. Please search for part number FC099 to order this attachment.

- diam: 6mm
- length: 164mm
- temp: 0 to 50°C
- max pressure: 0.1 bar

### Recommended use

- meat

### Ordering Information

pH SMART Electrode | FC232D  
supplied with BNC connector and 1m cable

#### Electrode Storage Solutions

electrode storage solution (500 mL) HI-70300L  
electrode storage solution (230 mL) HI-70300M

### Specifications

Product Code	FC232D
Reference	single, Ag/AgCl
Junction	open
Connector	BNC
Electrolyte	viscolene
Max Pressure	0.1 bar
Range	pH: 0 to 12; 0 to 50°C
Tip Shape	conic (6 x 10 mm)
Temperature Sensor	yes
Amplifier	yes
Body Material	PVDF
Recommended Use	meat

### DIN Connection



Authorised distributor

**Impex Produkter AS**  
Gamle Drammensvei 107  
1363 Høvik  
www.impex.no

info@impex.no  
Tel.: 22 32 77 20